

Client Update

August 2022

Page 1 of 2

Test Name	Test Code	Effective Date
-----------	-----------	----------------

IGF Binding Protein-3 **1316** **Immediately**

Specimen requirements for this test have been updated. Please see below for details.

Previous Test Information		Updated Test Information
Primary Container	SST	ALQS, ALQRD
Transportation Temp	Refrigerate	Strict Frozen

PTH Intact **0598** **8/8/2022**

REFERENCE RANGE UPDATE FOR ELMWOOD PARK, NJ LABORATORY

Due to a change in instrumentation, the Reference Range for test 0598 has been modified. See below for details.

Previous Test Information	Updated Test Information
Reference Range	18.4-80.1 pg/mL
	15.0-65.0 pg/mL

C-Telopeptide **T421** **8/8/2022**

Due to a change in instrumentation, the Reference Range for test T421 has been modified. See below for details.

Previous Test Information	Updated Test Information
Reference Range	Males(yrs) Range(pg/mL)
	<30 Not Established
	30-50 <585
	51-70 <705
	>70 <855
	Females Range(pg/mL)
	Premenopausal <574
	Postmenopausal 49-738
	Males(yrs) Range(pg/mL)
	<30 238-1019
	30-39 225-936
	40-49 182-801
	50-59 161-737
	60-69 132-752
	>69 118-776
	Females Range(pg/mL)
	<30 148-967
	30-39 150-635
	40-49 131-670
	50-59 183-1060
	60-69 171-970
	>69 152-858
	Premenopausal 136-689
	Postmenopausal 177-1015

Cord blood/Infants **TL74** **8/15/2022**

This is a new test code specific for cord blood samples and infants less than 4 months of age. This test code consists of ABO and Rh forward typing of infant red blood cells. Reverse typing is not recommended in this age group due to the immaturity of their immune system and the possibility of passively acquired anti-A or anti-B antibodies from the mother.

